Cerapulus S[®]

DIRECTIONS FOR USE

1. Spruing and waxing techniques

It is advisable to use sprues: 3-3.5mm in diameter.
Distribution channels: 4-5.0mm in diameter.
To make particularly thick castings or multi-element (4 or more units) bridges, it may be necessary to utilize larger sprue channels and/or to use sprue pins with spherical reservoir (5-7mm in diameter) about 2mm from the casting. This procedure or appropriate venting, dictated by personal experience, will reduce the risk shrinkage or lack of penetration during solidification.

2. Investing

Use only graphite and plaster-free phosphate-bonded investment compounds for high temperature applications. Carefully follow manufacturer's instructions. To permit free expansion of the compound during setting, line the inside of the ring with an adequate thickness of spacer to suit the ring size. Alternatively use a ring of your own choice which can expand freely to ensure accurate fitting castings.

3. Investment preheating

A guide line for preheating is to raise the temperature of the furnace slowly (5-7°C/minute) to about 300°C. Depending on the size and the number of rings, this temperature is held for 30-40 minutes, after which it is raised to 575°C at (5-7°C/minute) and held for 30-40 minutes. Finally it is raised to 950°C and again held for 30-40 minutes. Always refer to the manufacturer's instructions for the specific use of the casting investment.

4. CERAPULUS S: re-casting and over-heating

CERAPULUS S contains small amounts of certain chemical elements which favor the formation of special oxides that control oxidation kinetics and metal-ceramic bonding during porcelain firing. Alloy re-use or prolonged overheating during melting causes loss of these chemical elements by evaporation because of their low vapor pressure. This situation may have an adverse effect on the metal porcelain bond. Therefore for very precise and very thin castings it is advisable to use new crucibles and only new alloy. In the case of acrylic veneering, approximately 50% of previously used alloy can be employed. The alloy to be re-used must be carefully sandblasted (or preferably pickled with appropriate chemical agents) and cleaned by ultrasonic means or by a steam jet. The aim of these precautions is to ensure that the casting is not contaminated with impurities coming from the investment compound, or with sand-blasting residues or oxides which could affect the quality of the finished product.

5. Alloy melting and casting procedures

CERAPULUS S can be cast using the following methods:

- Oxygen-propane or oxygen butane flame.
 - High frequency induction furnaces working in air under vacuum or in a controlled atmosphere.
 - High-performance resistance casting furnaces capable of attaining at least 1450°C.
 - Centrifugal-casting or die-casting can be employed.
- It is not advisable to use either an oxy-acetylene flame or an electric arc as these tend to lead to overheating and carbon contamination of the molten alloy.

5.1 Flame casting

Preheat the empty crucible and adjust the gas flow to give a neutral flame, according to the instructions. The blue cone of the flame should be at least 3-5mm long. Set the torch about 40mm above the top of the crucible and rotate the flame continuously to obtain a uniform molten mass. Avoid breaking the oxide skin, since this may result in excessive oxidation. Add a small charge of flux. As soon as fissures begin to appear in the oxide skin start centrifugal-casting. It is recommended that the centrifuge arm be loaded with one revolution more than in the case of precious alloys to compensate for the lower density of CERAPULUS S. AVOID OVERHEATING THE ALLOY, which causes evaporation of the elements incorporated to facilitate ceramic bonding, results in spum splashing and leads to excessive oxidation which could greatly increase surface roughness of the casting.

5.2 High-frequency induction furnace casting

Set the temperature of the casting machine at about 1450°C (namely 130°C higher than the melting point of the alloy).

- Preheat the empty crucible.
- Place two or three metal ingots (10-15gr) in the crucible. Fix the investment ring preheated to 950°C.
- Start melting and add such metal as necessary, plus a small charge of flux. As soon as the oxide skin shows signs of opening up, start centrifugal or pressure casting.
- Leave the ring to cool to room temperature.

6. Sand-blasting and preliminary preparation

After casting is complete leave the ring to bench-cool slowly.

Remove the investment by sand-blasting with coarse quartz or aluminum (120-200µm) ensuring that the margins of the crowns are not sand-blasted directly.

Prepare the surfaces to be covered with ceramic abrasive stones or tungsten carbide burs and finish with a diamond bur applying the pressure in one direction only.

For acrylic and composite veneering, attach mechanical retention or use a silane process or other methods recommended by the manufacturer.

7. Preparation for porcelain veneering

For porcelain veneering, the alloy surface should be prepared as advised in point 6, and must be carefully sand-blasted again.

Prior to the final sand-blasting it is preferable to pickle the piece in a mixture of 10% HNO₃ and HF acid and 90% water for 30-40 minutes at room temperature, followed by careful rinsing in running water and/or by means of a steam jet. This precautionary operation removes the thin, adherent layer of oxide which is masked but not removed by the first coarse sand-blasting and also the residues of the chromium layer (1-5µm) which forms beneath the oxide film, and may not have been removed by the mechanical operations.

Sand-blasting must be done using pure new, un-recycled aluminum oxide in the 80 to 150µm size range at a relatively high pressure and for a sufficient length of time.

Never use materials that are excessively fine (around 50µm) because the grains tend to remain on the metal, thus interrupting the continuity of the porcelain-metal bonding. Clean the surface ultrasonically in distilled water and/or by means of a steam jet. Finally dip it for a few minutes in ethyl alcohol or ethyl acetate, possible in an ultrasonic tank.

From this moment the surface should be perfectly clean, with the appropriate surface whose chemical composition is identical with that of the super alloy, ready to receive the first layer of opaque or bonding agent under the best conditions to ensure an excellent porcelain-metal bond. Henceforth, the surface must be handled only with gloves or clean tongs.

8. Porcelain veneering

Do not plate with gold or other precious metals.

Do not perform any pre-oxidation treatment, typical of Au and Pd alloys or alloys containing beryllium.

With CERAPULUS S the opaque is applied directly in two layers to clean surface (see Point 7). The first, which is very dilute, is applied employing the wash brand technique and is then fired using exactly the temperature and conditions recommended by the porcelain manufacturer. The second layer of opaque must be thick enough to cover.

It is strongly recommended that opaque's with a high concentration of metal oxides be utilized, such as opaque pastes, since these favor development of the bond with beryllium-free non precious alloys. Alternatively, a suitable colour conditioner or bonding agent can be used instead of the first layer of opaque.

Cool the porcelain veneered case slowly to room temperature, following the scheduled cooling stages and the characteristics of the opaque employed.

Before every successive firing and cooling cycle, brush the case under running water and then re-clean in hot distilled water or with a steam jet.

CERAPULUS S is compatible with all kinds of porcelains available on the market.

9. Cleaning and Polishing

Clean and polish the porcelain veneered piece first by means of hard rubber wheels designed to remove scratches and the like from Ni-Cr and Co-Cr alloys. For final shiny polishing use water-based diamond paste (>5µm) and suitable brushes and felt pads.

Use an ultrasonic bath and/or a steam jet for final cleaning.

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